Shipping Biological - Blood Samples



See Product Number STP320

As recommend by Brock's Research and Academic Safety Officer, biological and blood samples should be shipped in similar containers as shown above. One supplier of this style of box is STP320.

http://www.saf-t-supply.com/shopper_catalog.asp?ID=2&ShopperRegionID=

These shipping containers are made for dry ice and can hold 4 x 50 ml test tubes: provided in kit. These are required packaging for shipping CL2 organisms which include human cell culture, blood samples, tissue that may contain a zoonotic (bat neural tissue). Staff members who are moving these materials themselves can reuse the containers or pick up the more durable re-useable version (often mailed back and forth).

Fedex offers some Pointers on Shipping: Clinical Samples, Biological Substance Category B(UN 3373) and Environmental Test Samples. To view them check out this link http://images.fedex.com/downloads/shared/packagingtips/pointers.pdf