

| Chemical Name |
|--|
| 2-Propanol - 4L Plastic Bottle |
| 30% Hydrogen Peroxide - 500ml |
| ACETONE -4L reagent |
| Acetone 20L |
| Acetone-D6 - ampule |
| ACETONE-D6 25G |
| Acetonitrile - HPLC Grade 4L |
| Acrylamide - Electrophoresis Grade - 500g |
| Agar - Bacteriological Grade 500g |
| Agarose - Biotechnology Grade - 100g |
| Ammonium Hydroxide - 2.5L |
| Benzene - Reagent Grade 4L |
| Boric Acid - Reagent grade - 500g |
| BRADFORD REAGENT - 500ML |
| Chloroform - HPLC grade - 4L |
| Chloroform-D3 - 100g |
| Chloroform-D3 Ampule |
| Citric Acid - reagent grade - 2.5kg |
| D2O - 100ml |
| D2O Deuterium oxide -ampule |
| DCM-D2 - 25G |
| Dextrose - reagent grade - 2.5kg |
| Dichloromethane - HPLC Grade - 4L |
| Dichloromethane - Reagent Grade - 20L |
| Diethyl Ether - 4L |
| Dimethyl Formamide - Reagent Grade - 4L |
| DIMETHYL SULFOXIDE 1L |
| DIMETHYL SULFOXIDE 4L |
| DMSO-D6 - 0.75ml ampule |
| DMSO-D6 - 50ml |
| Drierite - Blue - 1lb |
| Drierite - White - 1lb |
| ETHYL ACETATE - Reagent Grade - 20L |
| ETHYL ACETATE - Reagent Grade - 4L |
| Formaldehyde - 37% - 1L |
| Glacial Acetic Acid - reagent grade - 2.5L |
| Glycerol - reagent grade - 1L |
| Glycine - Biotechnology Grade - 500g |
| Hepes - reagent grade 100g |
| Hexanes - HPLC Grade - 4L |
| Hexanes - reagent grade - 20L |

| Chemical Name |
|---|
| Methanol - HPLC Grade - 4L |
| Methanol - 4L Plastic bottle |
| Methanol-D4 - 1ml ampule |
| Nitric Acid - 2.5L |
| Paraffin - Extra Heavy - 4L |
| Pentane - Reagent Grade - 4L |
| Peptone - 500g |
| Petroleum Ether - Reagent Grade - 4L |
| pH 2.00 Buffer - Clear - 500ml |
| pH 4.00 Buffer - Red - 500ml |
| pH 7.00 Buffer - Yellow - 500ml |
| pH 10.00 Buffer - Blue - 500ml |
| Phosphoric Acid - 2.5L |
| Polyvinyl Alcohol - Reagent Grade - 1kg |
| Potassium Carbonate |
| Potassium Chloride - Reagent - 2.5kg |
| Potassium Chloride - Reagent - 500g |
| Potassium Hydroxide - 2.5kg |
| Propionic Acid - 2.5L |
| Soda Lime - Indicating 500g |
| Sodium Acetate Anhydrous - 500g |
| Sodium Bicarbonate - 2.5kg |
| Sodium Bicarbonate - 500G |
| Sodium Carbonate - 500g |
| Sodium Chloride - 2.5kg |
| Sodium Citrate - Reagent - 1kg |
| Sodium Hydroxide Pellets - 2.5kg |
| sodium Phosphate Dibasic |
| sodium phosphate monobasic - 500g |
| Sodium Sulfit - 500g |
| Sodium Sulphate - 2.5kg |
| Sterile Water - 1000ml |
| Sucrose - 2.5kg |
| Sulphuric Acid - 2.5L |
| Tetrahydrofuran - Reagent Grade - 4L |
| Toluene - Reagent Grade - 4L |
| Tris(hydroxymethyl)aminomethane - 1kg |
| Tryptone - Bacteriological Grade 500g |
| Water - HPLC Grade - 4L |
| XYLENES |
| Yeast Extract |

| |
|----------------------------|
| Hydrochloric Acid - 2.5L |
| Hydrochloric Acid - 500ml |
| Iodine - 100g |
| Magnesium Sulphate - 2.5kg |

